



# BIKO SPORTS BANDS

## SIZING AND FITTING YOUR DOG FOR THE BIKO SPORTS BANDS

Accurate measurements, whether in inches or centimeters, ensure a better fit.

Biko Sports Bands are available in sizes 3 through 6. Sizes 1 and 2 for smaller breeds, and size 7 for Giant breeds can be made upon request at no additional costs.

BAND SIZING CHART		
Band Size	Centimeters	Inches
1	40 to 45	15.75 to 17.75
2	46 to 51	18 to 20
3	52 to 59	20.5 to 23.25
4	60 to 67	23.5 to 26.5
5	68 to 75	26.75 to 29.5
6	76 to 82	30 to 32.25
7	83 to 90	32.5 to 35.5

The padded foot cuffs are available in sizes 1 through 4.

FOOT CUFF SIZING CHART		
Cuff Size	Centimeters	Inches
1	5 to 7	2 to 2.75
2	7.5 to 10	3 to 4
3	10.5 to 13	4.25 to 5
4	13.5 to 18	5.25 to 7

### THE PROCEDURE FOR FITTING THE SPORTS BANDS IS SIMPLE AND STRAIGHTFORWARD:

- Step 1:** Put the Roman harness on the dog, D ring lying flat.
- Step 2:** Place the padded cuffs (with Velcro fasteners) around the metatarsus on each hind limb, just below the hock.
- Step 3:** Clip the Biko PR Bands (the end without extenders) to the D ring of the harness.
- Step 4:** Clip the end with the extenders to the foot cuffs.



Typical breeds for size 1: Beagle and Pug



Typical breeds for size 2: Miniature Schnauzer and French Bulldog



Typical breeds for size 3: Shetland Sheepdog and Pit Bull Terrier



Typical breeds for size 4: Border Collie and Vizsla



Typical breeds for size 5: German Shepherd and Labrador Retriever



Typical breeds for size 6: Bernese Mountain Dog and English Mastiff



Typical breeds for size 7: Newfoundland and Great Dane



The animal mobility & fitness company

We can easily accommodate special orders placed for dogs under 40 centimeters (15.75 inches) or over 90 centimeters (36 inches).



[www.aniMotionProducts.com](http://www.aniMotionProducts.com)